


EASY UPGRADE SOLUTION.

FluidFM® ADD-ON FOR YOUR EXISTING AFM INSTRUMENT



EASY UPGRADE SOLUTION.

Example of a FluidFM ADD-ON installed on a Nanosurf FlexAFM.

sub-nm

XYZ-RESOLUTION

pN

FORCE CONTROL

100% AFM


CAPABILITIES

100% FluidFM®

APPLICATIONS

Empower your AFM beyond imaging with FluidFM. The FluidFM ADD-ON solutions from Cytosurge allow you to enter the world of FluidFM using your existing atomic force microscope.

The innovative upgrade solution opens the door to a wide range of new applications made possible by FluidFM technology. From single cell biology to material sciences – harness the possibilities provided by FluidFM using the AFM you already own.



CLIP IT ON.

Close-up of the Nanosurf FlexAFM scan head with a mounted FluidFM probe.

MULTIPLY YOUR POSSIBILITIES WITH THE FluidFM® ADD-ON.

FluidFM technology extends the application scope of classical AFM way beyond imaging and force spectroscopy. Applications ranging from single cell isolation to colloidal spectroscopy measurements become readily possible with this easy upgrade solution.

FluidFM technology enables your research at the nanoscale to become faster, more powerful and economical. Discover a multitude of proven FluidFM applications ranging from biophysics to material sciences and beyond. The possibilities provided by FluidFM are endless. Unleash the true power of your AFM with the FluidFM ADD-ON from Cytosurge.

YOUR REQUIREMENTS. OUR SOLUTIONS.

Developed by scientists for scientists; we know that your research is as unique as yourself. Because we want to offer you the best possible solutions for your specific requirements we have created a bespoke selection of packages that allows you to choose the perfect product to match your needs.

FluidFM® ADD-ON Silver	FluidFM® ADD-ON Gold	FluidFM® ADD-ON Platinum
---------------------------	-------------------------	-----------------------------

The **Silver** package is the ideal choice for FluidFM experts looking for an economical start into the world of FluidFM. The Silver package offers 12 months warranty and it includes the Silver FluidFM Starter kit containing 3 probes.


The **Gold** package is our most popular choice, maximizing your value while playing nicely on your budget. The Gold package includes everything to ensure a smooth start into the world of FluidFM applications. You can benefit from 24 months warranty. The Gold package includes the Gold FluidFM Exploration kit containing a selection of 10 micropipettes and nanopipettes of various stiffnesses and apertures.

CONTACT US.

We and our distributors offer comprehensive support for FluidFM® customers. Please visit the Cytosurge Help Center in order to access the FluidFM® user community: support.fluidfm.com

CYTOSURGE®

CYTOSURGE AG, SÄGEREISTRASSE 25, 8152 GLATTBRUGG, SWITZERLAND
PHONE +41 43 544 87 00, FAX +41 43 544 87 09, WWW.CYTOSURGE.COM



EMPOWER YOUR AFM.

Standard AFM versus FluidFM technology.

Our **Platinum** package is the best choice for you if you want to experience the full power of the FluidFM ADD-ON solution. The integrated video camera* allows you to stay in full control of your experiment at all times. Thanks to the smart probe filling assistant**, getting your experimental setup ready has never been easier. With the Platinum package, you have 36 months warranty and you can take advantage of the comprehensive Platinum FluidFM Discovery kit or the FluidFM Explorer kit and additional 10 probes of the same type of your choice. These two options enable you to select the probes kit that best fits your needs.

The FluidFM ADD-ON for AFM is compatible with scan heads from major AFM manufacturers. Please contact us to inquire if your AFM is supported.

* Requires inverted optical microscope with free access port. Please inquire with us to check if your setup is compatible.

** Software tool for real-time, reliable detection of probe filling.

